

Agrobacterium Transformation and Competent Cell Preparation

Monday, January 07, 2013

3:59 PM

Sean Weise

Agrobacterium Transformation

Materials:

Gene-Pulse Cuvettes, 0.2cm (BIO-RAD #1652086)

LB

Spectinomycin

Rifampicin

LB plates with antibiotic

Protocol:

1. Dilute plasmid to 15 ng/ μ L (if the [salt] is too high, the cells will be killed by the electrical impulse)
2. Mix 2 μ L diluted plasmid to 50 μ L GV3101 agrobacterium in a 1.5mL tube
3. Transfer mixture into a cuvette that has been chilled for 20 min @ -20°C by adding to the side & then tapping down to eliminate air bubbles
4. Set the electroporator to 2.00 V (leave all other settings at default: 200 Ω , capacitance extender 250 μ FD, capacitance 25 μ FD)
5. Snap cuvette into proper orientation, push in & hold both pulse buttons until beep (may take a few seconds) the time constant should read about 4.7
6. Add 1mL LB or SOC to cuvette & pipette up & down to mix, transfer 1 ml microfuge tube
7. Incubate at 28°C (cells will die at 37°C) for 2-3 hours before plating ~70 μ L on LB+spectinomycin+rifampicin plates
8. Wrap plates with micropore & incubate at 28°C for 2-3 days, it really takes at least 2 days, even at 28°C so don't be disappointed if you don't see colonies the next day, be patient
9. Make liquid cultures from colonies for PCR testing

Agrobacterium Competent Cell Preparation

Materials

LB plates with 30 μ g/ml Gen (2 plates is enough)

2 1 L culture flasks

2 sterile 500 ml centrifuge bottles

1.5 - 2 L sterile dH₂O (should be cold)

50 ml sterile 10% glycerol in dH₂O (should be cold)

Protocol:

Note: The Agrobacterium we use is GV3101, which has both resistance to Rif (10 μ g/ml, in genome), Gen (30 μ g/ml, in helper plasmid) This protocol is often done with any antibiotics but you can add Gen to plates and starter liquid cultures to be one the safe side.

1. From a frozen glycerol stock grow a 20 ml culture of Agrobacterium overnight at room temperature
2. The following morning dilute the culture 64,000 fold by adding 0.5 μ l of culture to 2 ml of H₂O making a **4000 fold** dilution. From the 4000 fold dilution take 0.5 μ l and add to 2 ml H₂O making a **16,000 fold** dilution, from the 16,000 fold dilution take 0.5 μ l and add to 2 ml H₂O making a **64,000 fold** dilution
Note: 0.5 μ l of culture into 50 ml H₂O makes a 100,000 fold dilution and you can use 10 μ l of this on a plate
3. Add 5 μ l of the 64,000 fold dilution to an LB plate and thoroughly spread the culture on the plate, incubate at room temp or 28°C for 1 - 3 days
4. Once you have single colonies on a plate, pick a colony and start a 20 ml overnight culture in LB with Gen at 28°C , you can start this overnight culture early in the day as Agro grows slow and you want a thick culture for the next day. Make sure to take LB out of the fridge so that it is at room temperature for the 500 ml cultures the next day
5. Inoculate 2 500 ml flasks (use 1 liter flasks) of LB with 9 ml of the overnight culture, do not use antibiotic here, Make sure to start these cultures early in the day (if possible < 8:00am) as Agro grows very slowly, a culture started around

- 8:00am will be ready hopefully by 3:00pm
6. Grow cells to an OD_{600 nm} of 0.5 - 1
 7. Precool Centrifuge with 500 ml bottle adaptors to 4°C
 8. When cells are ready to harvest chill flasks on ice for 15 - 30 minutes
 9. Centrifuge at 4000 g in bottles for 15 minutes at 4°C
 10. Remove as much of the supernatant as possible and resuspended the pellet in each bottle in 500 ml of ice cold water, use graduations on side of centrifuge bottle to measure volume
 11. Centrifuge at 4000 g for 15 min
 12. Remove supernatant and resuspended the pellet in each bottle in 250 ml of ice cold water
 13. Centrifuge at 4000 g for 15 min
 14. Remove supernatant and resuspended the pellet in each bottle in 10 ml ice cold 10% glycerol
 15. Combine contents of bottle into 1 50 ml Falcon tube (It will stand up the centrifugation force)
 16. Centrifuge at 4000 g for 15 min
 17. Remove supernatant and resuspend the pellet in 2 - 3 ml ice cold 10% glycerol. The cell concentration should be about $1 - 3 \times 10^{10}$ cells/ml
 18. Store in 1.5 ml microfuge tubes in 50 µl aliquots, freeze in LN₂, store in -80°C, cells are good for at last 6 months

Adapted from protocol by Sanjay

Agrobacterium Transformation Reagents

LB Liquid, 200 ml

In 200 ml dH₂O
 + 5 g LB Broth (Difco, Luria-Bertani)
 Autoclave for 20 min

LB Plates, 200 ml

In 200 ml dH₂O
 + 5 g LB Broth (Difco, Luria-Bertani)
 + 3 g Agar (Fishger, BP 1423-500)
 Autoclave for 20 min, swirl immediately, allow to cool to touch
 Pour plates (≈ 35 ml) label, store at 4°C

4 mg/ml Ampicillin, 10 ml

In 10 ml of sterile dH₂O
 + 0.04g Ampicillin (Roche 70043429)
 Filter sterilize using 0.45 µm syringe filter (Millipore Millex HA SLHA033SS)
 Aliquot into sterile 1.5 ml microfuge tubes and store at -20°C

10 mg/ml Spectinomycin, 10 ml

In 10 ml sterile dH₂O
 + 0.1 g Spectinomycin
 Filter sterilize using 0.45 µm syringe filter (Millipore Millex HA SLHA033SS)
 Aliquot into sterile 1.5 ml microfuge tubes and store at -20°C

34 mg/ml Rifampicin, 10 ml

In 10 ml sterile dH₂O
 + 0.34 g Rifampicin
 Filter sterilize using 0.45 µm syringe filter (Millipore Millex HA SLHA033SS)
 Aliquot into sterile 1.5 ml microfuge tubes and store at -20°C in the dark, rifampicin is light sensitive

Antibiotic	Stock (mg/ml)	Working (µg/ml)	µl of Stock in 5 ml	10 ml	20 ml	200 ml
Ampicillin	10	50	12.5	25	50	500
Spectinomycin	10	100	50	100	200	2,000
Rifampicin*	34	150	22	44	88	882
Gentamicin	10	30	15	30	60	600

* light sensitive: store stock solutions and plates in the dark

3:59 PM

10% Glycerol, 200 ml

In 180 ml dH₂O

+ 20 ml glycerol

Autoclave for 20 - 30 min, cool and store in cold room