

# Bacterial Transformation and Plasmid Prep

Monday, November 07, 2011

11:15 AM

Sean Weise

## Growing up plasmid in *E. coli* for glycerol stocks and plasmid stock

### Materials:

DH5 $\alpha$  Competent *E. coli* strain, -80°C

Precision H<sub>2</sub>O bath at 42°C

Sterile 2 ml microfuge tubes

Plasmid (5 ng /  $\mu$ l)

LB plates with antibiotic

LB liquid

sterile cryoviles

sterile 80% glycerol

LN<sub>2</sub>

QIAprep Spin Miniprep Kit (Quiagen, 27104)

P1 buffer from above kit, 4°C

1.5 ml microfuge tubes sterile, with top cut off

0.5 ml sterile microfuge tubes

### Methods:

#### Transforming and selecting *E. coli* with your plasmid

1. Thaw DH5 $\alpha$  *E. coli* cells on ice, must be kept on ice cannot come to room temp!
2. Turn on Precision H<sub>2</sub>O bath to reach 42°C, put LB in bath to warm up
3. In a 1.5 ml microfuge tube, mix 50  $\mu$ l competent cells with 2 - 4  $\mu$ l plasmid (10 ng DNA), treat gently don't vortex
4. Incubate on ice for 0.5 hr
5. Heat shock by putting in floating rack in water bath at 42°C for 20 s
6. Immediately put back on ice for 2 minutes
7. Add 450  $\mu$ l of prewarmed LB (if SOC is available use SOC)
8. Incubate tubes for 1 hour at 37°C, 225 rpm on shaker (can be slower)
9. Put  $\approx$  80 - 150  $\mu$ l of *E. coli* on LB agar plates with antibiotic  
You can save the extra *E. coli* at 4°C overnight in case your plates don't work
10. Incubate the plates over night at 37°C

#### Growing up lots of *E. coli* with plasmid

1. Prepare a small 50 ml sterile glass flask with 5 ml of LB and antibiotic (same concentration as on plates)
2. Take 1 colony from plate with 200  $\mu$ l pipette tip and place in flask  
Grow up 2 cultures from 2 separate colonies per plate, in case one doesn't work
3. Incubate overnight at 37°C, 200 rpm on shaker

#### Glycerol Stock

1. To make glycerol stock prepare a sterile cryovile by adding 800  $\mu$ l of 80% sterile glycerol
2. To the cryovile add 1 ml of *E. coli* liquid culture (try to catch culture in exponential phase of growth)
3. Snap freeze in LN<sub>2</sub>

#### Mini Prep to isolate plasmid

##### If using Miniprep kit for first time (i.e. a new miniprep box)

1. Add the RNase A solution to buffer P1
  2. Add lysis blue to buffer P1
  3. Check box on buffer P1 that RNase A has been added, and label bottle top that lysis blue has been added
  4. Place date on bottle (**Only good for 6 months**)
  5. Add ethanol to buffer PE, will say on bottle how much ethanol to add
  6. Label buffer PE that ethanol has been added
- 
1. Place 2 ml of *E. coli* liquid culture into a 2 ml microfuge Tube and centrifuge for 10 min
  2. Discard supernatant in biohazard waste and resuspend pellet in 250  $\mu$ l P1 buffer (in 4°C) break up pellet with pipette tip  
\*\* P1 buffer is good for 6 months so watch date on bottle \*\*
  3. Add 250  $\mu$ l of P2 cell lysis buffer, mix by inverting until liquid is all blue  
\*\* Do not allow to remain in P2 buffer for more than 5 minutes \*\*
  4. Add 350  $\mu$ l of buffer N3, mix by inverting until liquid is clear and SDS precipitate forms
  5. Centrifuge at max speed for 10 minutes, white pellet will form this is genomic DNA and junk

6. Pipette supernatant onto spin column and discard tube with pellet in biohazard waste
7. Centrifuge for at max speed for 1 minute, plasmid will stick to membrane  
 Don't need to wash is PB buffer after this step when using DH5 $\alpha$  E. coli  
 If needed add 0.5 ml PB buffer centrifuge 1 min discard flow through  
 - This removes trace nuclease activity with using EnAt strains
8. Discard flow through in biohazard waste and add 750  $\mu$ l of PE buffer  
 \*\* Make sure ethanol has been added to the PE buffer \*\*
9. Centrifuge for 1 minute and discard flow through
10. Centrifuge again for 1 minute to remove residual liquid and discard flow through
11. Transfer spin column to a new sterile 1.5 ml tube with top cut off
12. Add 50  $\mu$ l EB or DNase RNase free H<sub>2</sub>O to center of membrane, let stand for 1 minute
13. Centrifuge for one minute, this is your purified plasmid

**Check plasmid quality, concentration on spectrophotometer**

1. Use nano-drop spectrophotometer, choose nucleic acids
2. Blank with EB buffer
3. 230 nm = polyphenolics, 260 = DNA, 280 = protein  
 $260/280 > 1.8$   
 $230/260 > 2$

**Reagents for Growing E.Coli:**

**LB Plates, 200 ml**

In 200 ml dH<sub>2</sub>O  
 + 5 g LB Broth (Difco, Luria-Bertani)  
 + 3 g Agar (Fishger, BP 1423-500)  
 Autoclave for 20 min, swirl immediately, allow to cool to touch  
 Pour plates ( $\approx$  35 ml) label, store at 4°C

**LB Liquid, 200 ml**

In 200 ml dH<sub>2</sub>O  
 + 5 g LB Broth (Difco, Luria-Bertani)  
 Autoclave for 20 min

**80% Glycerol, 80 ml**

In 16 ml dH<sub>2</sub>O  
 + 64 ml Glycerol  
 Autoclave for 20 min

| Antibiotic                 | Stock (mg/ml) | Working ( $\mu$ g/ml) | $\mu$ l of stock in 5 ml | $\mu$ l of stock in 20 ml | $\mu$ l of stock in 40 ml | ml of stock in 200 ml | ml of stock in 500 ml | ml of stock in 1 L |
|----------------------------|---------------|-----------------------|--------------------------|---------------------------|---------------------------|-----------------------|-----------------------|--------------------|
| Ampicillin                 | 10            | 50                    | 25                       | 100                       | 200                       | 1                     | 2.5                   | 5                  |
| Chloramphenicol (methanol) | 10            | 20                    | 10                       | 40                        | 80                        | 0.4                   | 1                     | 2                  |
| Kanamycin                  | 10            | 25                    | 12.5                     | 50                        | 100                       | 0.5                   | 1.25                  | 2.5                |
| Rifampicin* (methanol)     | 30            | 150                   | 25                       | 100                       | 200                       | 1                     | 2.5                   | 5                  |
| Spectinomycin              | 10            | 100                   | 50                       | 200                       | 400                       | 2                     | 5                     | 10                 |
| Streptomycin               | 10            | 30                    | 15                       | 60                        | 120                       | 0.6                   | 1.5                   | 3                  |
|                            |               |                       |                          |                           |                           |                       |                       |                    |

Note: If using multiple antibiotics you can cut the concentration of each by half

\* Rifampicin is light sensitive cover with aluminum foil

**10 mg/ml Antibiotic, 10 ml**

In 10 ml of sterile dH<sub>2</sub>O  
 + 0.1g Antibiotic  
 Filter sterilize using 0.2 - 0.45  $\mu$ m syringe filter  
 Aliquot into sterile 1.5 ml microfuge tubes and store at -20°C