

DNase Treatment using RQ1 (Promega)

Wednesday, November 02, 2011
11:57 AM



RQ1 RNAs...

I. DNase Treatment of RNA Samples Prior to RT-PCR

1. Set up the DNase digestion reaction as follows:

RNA in water or TE buffer 1–8 μ l

RQ1 RNase-Free DNase 10X Reaction Buffer 1 μ l

RQ1 RNase-Free DNase (see Note 1) 1u/ μ g RNA

Nuclease-free water to a final volume of 10 μ l

2. Incubate at 37°C for 30 minutes.

Note: If analyzing RNA samples by gel electrophoresis, perform a phenol:chloroform extraction and ethanol precipitation before loading the samples on the gel, because salts in the RQ1 DNase Reaction Buffer and Stop Solution cause aberrant migration of RNA on gels. Steps 3 and 4 may be omitted if a phenol:chloroform extraction is performed.

3. Add 1 μ l of RQ1 DNase Stop Solution to terminate the reaction.

4. Incubate at 65°C for 10 minutes to inactivate the DNase.

5. Add all, or a portion of, the treated RNA to the RT-PCR reaction. See the Access RT-PCR System(a) Technical Bulletin #TB220.

Notes:

1. Use 1 unit of RQ1 RNase-Free DNase per microgram of RNA. For smaller amounts of RNA, use 1 unit of RQ1 RNase-Free DNase per reaction.

2. The RQ1 RNase-Free DNase digestion contains a final concentration of 10mM MgSO₄. When adding the DNase-treated RNA to an RT-PCR reaction, carryover of magnesium must be considered. For example, the addition of 1 μ l of treated RNA to a 50 μ l RT-PCR reaction will raise the magnesium concentration by 0.2mM, and the addition of 5 μ l of treated RNA will raise the magnesium concentration by 1mM. The requirement for magnesium may be different in the RQ1 DNase digestion step and in the amplification reaction.

- RQ1 DNase activity increases as the Mg²⁺ concentration increases up to 5–10mM. At a concentration of 1mM Mg²⁺, RQ1 DNase is expected to be at least fourfold less active than at the optimal Mg²⁺ concentration.

- For some templates, the yield from the amplification reaction is highly dependent on the Mg²⁺ concentration, and the optimal Mg²⁺ concentration may be as low as 1mM.

If an increased Mg²⁺ concentration is not tolerable in the amplification reaction, the following alternatives may be used.

- The RQ1 RNase-Free DNase 10X Reaction Buffer may be diluted 1:10 with 400mM Tris (pH 8.0), 10mM CaCl₂ prior to DNase digestion. (Note that, under these conditions, the RQ1 DNase will be approximately fourfold less active than under standard reaction conditions.)

- An alternative DNase reaction buffer may be used (such as the RT or PCR

reaction buffer) if that buffer contains at least 1mM Mg²⁺.

- The RNA sample may be diluted in water prior to RT-PCR allowing dilution of the MgSO₄ to a concentration that is compatible with this application.
- The RNA may be purified with a standard phenol:chloroform extraction followed by an ethanol precipitation.