

Growth Media

Tuesday, March 27, 2012
2:20 PM

LB Liquid

200 mL dH₂O
+5g LB broth (Acumedia #7279)
Autoclave for 25 minutes

LB Plates

200 mL dH₂O
+5g LB broth (Acumedia #7279)
+3g Agarose (Invitrogen #16500-500)
Autoclave for 25 minutes and allow to cool to touch. Add appropriate amount of antibiotics and pour ~35mL/plate in a sterile workspace. Store at 4°C

MM3 Media for Batch Cultures

Adjust Glucose and yeast as needed
Seed culture we typically use 5 g/L glucose, 5g/L yeast extract
Experimental culture we typically use 5 g/L glucose, 1 g/L yeast extract

Components	1 L	250 ml	150 ml	
dH ₂ O	600 ml	150 ml	90 ml	
50% Glucose (5 g L ⁻¹ , 27.8 mM batch)	8.2 ml	2.05 ml	1.25 ml	
1 M MOPS	100 ml	25 ml	15 ml	
10X Medium concentrate	100 ml	25 ml	15 ml	
10X Phosphate buffer	100 ml	25 ml	15 ml	
20X CaCl ₂	50 ml	12.5 ml	7.5 ml	
100X Trace Elements	10 ml	2.5 ml	1.5 ml	
1 mg/ml Thiamine (2 mg L ⁻¹ batch)	2 ml	0.5 ml	0.3 ml	
100 g/L Yeast extract (1 g L ⁻¹ batch)	10 ml	2.5 ml	1.5 ml	

pH to 7.0 using NaOH or HCl (usually not necessary)
Adjust volume to final volume
Filter sterilize and store at 4°C

Concentrated Solutions

50% Glucose Solution 250 mL

In 250mL dH₂O
+ 125 g Glucose (fw 180) 5.76 M stock, 94.5 mM batch
Filter sterilize or autoclave

Yeast Extract 100 g/L

In 150 ml

+ 15 g of yeast extract
Filter sterilize

1 mg/ml Thimine Solution

In 50 ml H₂O
+ 0.05 g Thiamine HCl
Filter sterilize and store at 4°C

Medium Concentrate 10x

Dissolve following in 850mL dH₂O:
+ 50 g (NH₄)₂SO₄ (fw 132.1) 380 mM stock, 38 mM batch
+ 5 g MgSO₄ (fw 120.3) 41.5 mM stock, 4 mM batch
Adjust to 1L after dissolved for final volume

Phosphate Buffer 10x

Dissolve following in 1 L dH₂O:
+ 4 g KH₂PO₄ (fw 136.09) 29 mM stock, 2.9 mM batch
+10 g Na₂HPO₄ (fw 141.96) 71 mM stock, 7.1 mM batch
Adjust to pH of 7 with NaOH or HCl

Calcium chloride Solution 20x

In 1L dH₂O
+ 2 g CaCl₂ (fw 147) 13.5 mM stock, 0.7 mM batch

MOPS Buffer 1M

In 600 ml dH₂O
+ 209.3 g MOPS (fw 209.26) 1 M stock, 100 mM batch
Adjust to 1 L after adding MOPS for final volume
Adjust to pH 7.0 with NaOH or HCl

Trace Elements Solution 100x

NOTE: The weight for chemicals with different hydration levels should be adjusted on molecular weight base.

In 1L of dH₂O
+10 g Citric acid monohydrate (fw 210.14) 47.5 mM stock, 475 µM batch
+2 g MnSO₄*2H₂O (fw 183)* 10.9 mM stock, 109 µM batch
+0.5 g FeSO₄*7H₂O (fw 277.9) 1.8 mM stock, 18 µM batch
+0.2 g ZnSO₄*7H₂O (fw 280.54) 710 µM stock, 7.1 µM batch
+ 0.05 g CuSO₄*5H₂O (fw 249.6) 200 µM stock, 2 µM batch
+0.05 g Na₂MoO₄*2H₂O (fw 241.95) 207 µM stock, 2.1 µM batch
+ 0.05 g CoCl₂*6H₂O (fw 237.93) 210 µM stock, 2.1 µM batch
Filter sterilize and store at 4°C

MM1 for Batch Cultures

Materials:

Glucose
(NH₄)₂SO₄
MgSO₄
KH₂PO₄
Na₂PO₄

CaCl₂*2H₂O
Citric acid monohydrate
MnSO₄*2H₂O
FeSO₄*7H₂O
ZnSO₄*7H₂O
CuSO₄*5H₂O
NaMoO*2H₂O
Yeast extract
Thiamine
Filter sterilization flasks, 1 L and 150 ml

Stocks:

50% Glucose Solution, 250 mL

In 250mL dH₂O
+ 125 g Glucose (fw 180) 5.76 M stock, 94.5 mM batch
Filter sterilize or autoclave

10X Medium concentrate, 1L

Dissolve following in 850mL dH₂O:
+ 50 g (NH₄)₂SO₄ (fw 132.1) 380 mM stock, 38 mM batch
+ 10 g MgSO₄ (fw 120.3) 83 mM stock, 8 mM batch
Adjust to 1L after dissolved for final volume
Autoclave

10X Phosphate Buffer, 1 L

Dissolve following in 850mL dH₂O:
+ 50 g KH₂PO₄ (fw 136.09) 367 mM stock, 36.7 mM batch
+120 g Na₂HPO₄ (fw 140.95) 850 mM stock, 85 mM batch
Titrate to a pH of 7 with 5N NaOH or H₂SO₄
*NOTE: The Na₂HPO₄ takes a long time to dissolve, add to solution in 30 g increments
Autoclave

20X CaCl₂ Solution, 1 L

In 1L dH₂O
+ 4 g CaCl₂ (fw 147) 27 mM stock, 1.36 mM batch
Autoclave

MM1 Media for Batch Cultures Methods:

100X Trace Elements Solution, 1 L

NOTE: The weight for chemicals with different hydration levels should be adjusted on molecular weight base.

In 1L of dH₂O
+10 g Citric acid monohydrate (fw 210.14) 47.5 mM stock, 475 μM batch
+2 g MnSO₄*2H₂O (fw 183)* 10.9 mM stock, 109 μM batch
+0.5 g FeSO₄*7H₂O (fw 277.9) 1.8 mM stock, 18 μM batch
+0.2 g ZnSO₄*7H₂O (fw 280.54) 710 μM stock, 7.1 μM batch
+ 0.02 g CuSO₄*5H₂O (fw 249.6) stock, mM batch
+0.02 g Na₂MoO₄*2H₂O (fw 241.95) 82.6 μM stock, 0.826 μM batch
Autoclave

100 g/L Yeast Extract

In 150 ml

+ 15 g of yeast extract

Filter sterilize

1 mg/ml Thimine SolutionIn 50 ml H₂O

+ 0.05 g Thiamine HCl

Filter sterilize and store at 4°C

MM1 Media for Batch Cultures

Components	1 L	250 ml	150 ml	
dH ₂ O	600 ml	150 ml	90 ml	
50% Glucose (10 g L ⁻¹ , 55.5 mM batch)	16.4 ml	4.1 ml	2.5 ml	
10X Medium concentrate	100 ml	25 ml	15 ml	
10X Phosphate buffer	100 ml	25 ml	15 ml	
20X CaCl ₂	50 ml	12.5 ml	7.5 ml	
100X Trace Elements	10 ml	2.5 ml	1.5 ml	
1 mg/ml Thiamine (2 mg L ⁻¹ batch)	2 ml	0.5 ml	0.3 ml	
100 g/L Yeast extract (1 g L ⁻¹ batch)	10 ml	2.5 ml	1.5 ml	

pH to 7.0 using NaOH or HCl

Adjust volume to final volume

Filter sterilize and store at 4°C