Lowry Protein Assay

Tuesday, June 19, 2012 3:48 PM

METHOD 1 - 96 Well Plate

Materials

- Modified Lowry Protein Assay Kit (contents: Modified Lowry Protein Assay Reagent, Folin & Ciocalteu's Phenol Reagent, Albumin Standard, ampules) (From Kit: Thermo Scientific cat# 23240)
- 2 mL test tubes (10)
- Elution buffer
- Autoclaved water
- COSTAR (#3595)96 well plate and cover
- Multi-channel pipette
- Spectrophotometer

Methods

NOTE: Do not use Sarstedt Plates

1. Use one Albumin Standard (2.0 mg/mL) as stock and use Elution Buffer as Diluent Prepare the following BSA Standards Accordingly. Store in -20°.

Vial (name)	Volume of Diluent	Volume and Source of BSA	Final BSA Concentration
Α	125 μL	375 μL of stock	1,500 μg/mL
В	312.5 μL	312.5 μL of stock	1,000 μg/mL
С	155 μL	155 μL of vial A dilution	750 μg/mL
D	312.5 μL	312.5 μL of vial B dilution	500 μg/mL
Е	312.5 μL	312.5 μL of vial D dilution	250 μg/mL
F	312.5 μL	312.5 μL of vial E dilution	125 μg/mL
G	400 μL	100 μL of vial F dilution	25 μg/mL
Н	400 μL	100 μL of vial G dilution	5 μg/mL
I	400 μL	100 μL of vial H dilution	1 μg/mL
J	500 μL	0	0 μg/mL = blank

- 2. Prepare the Folin-Ciocalteu Reagent by diluting the 2X supplied reagent 1:1 with autoclaved water. Each test replicate requires 20 μ L of 1X Folin-Ciocalteu Reagent when used with a 96-well plate. The diluted reagent is unstable and must be made on the same day as use
- 3. Pipette 40 µL of each standard and unknown sample into a separate well
- 4. Add 200 μ L of Modified Lowry reagent to each well at nearly the same moment using a multichannel pipette
- 5. Immediately mix microplate with spectrophotometer for 30 seconds
- 6. Cover the plate with plastic lid and incubate at room temperature for exactly ten minutes
- 7. Add 20 μ L of 1X Folin-Ciocalteu Reagent to each well using a multi-channel pipette. Immediately mix with spectrophotometer for 30 seconds
- 8. Cover plate with plastic lid and incubate at room temperature for 30 minutes
- 9. Measure the absorbance at 750 nm on a plate reader. Subtract the average 750 nm absorbance value on the Blank standard replicates from the 750 nm value of the other individual standards and unknown sample replicates

10. Prepare a standard curve by plotting the average Blank-corrected 750 nm values for each BSA standard vs. its concentration in μ g/mL. Use the standard curve to determine the protein concentration for each unknown sample.

METHOD 2 – Cuvettes

Materials

- Modified Lowry Protein Assay Kit (contents: Modified Lowry Protein Assay Reagent, Folin & Ciocalteu's Phenol Reagent, Albumin Standard, ampules)
- 2 mL test tubes (10)
- Elution buffer
- Autoclaved water
- Polystyrene Cuvettes
- Spectrophotometer

Methods

1. Use two Albumin Standard (2.0 mg/mL) as stock and use Elution Buffer as Diluent Prepare the following BSA Standards Accordingly. Store in -20°.

Vial (name)	Volume of Diluent	Volume and Source of BSA	Final BSA Concentration
Α	250 μL	750 μL of stock	1,500 μg/mL
В	625 μL	625 μL of stock	1,000 μg/mL
С	310 μL	310 μL of vial A dilution	750 μg/mL
D	625 μL	625 μL of vial B dilution	500 μg/mL
E	625 μL	625 μL of vial D dilution	250 μg/mL
F	625 μL	625 μL of vial E dilution	125 μg/mL
G	800 μL	200 μL of vial F dilution	25 μg/mL
Н	800 μL	200 μL of vial G dilution	5 μg/mL
1	800 μL	200 μL of vial H dilution	1 μg/mL
J	1000 μL	0	0 μg/mL = blank

- 2. Prepare the Folin-Ciocalteu Reagent by diluting the 2X supplied reagent 1:1 with autoclaved water. Each test replicate requires 100 μ L of 1X Folin-Ciocalteu Reagent when used with cuvettes. The diluted reagent is unstable and must be made on the same day as use
- 3. Pipette 50 µL of each standard and unknown sample into a separate well
- 4. Add 1 mL of Modified Lowry reagent to each cuvette at 15-second intervals. Mix well and incubate at room temperature for exactly 10 minutes.
- 5. Maintain the 15-second interval between cuvettes by adding $100 \mu L$ of 1X Folin-Ciocalteu Reagent to each cuvette exactly at the end of each cuvette's 10-minute incubation period. Mix well.
- 6. Cover and incubate at room temperature for 30 minutes
- 7. Measure the absorbance at 750 nm with a spectrophotometer. Zero using a cuvette filled with only water. Subtract the average 750 nm absorbance value on the Blank standard replicates from the 750 nm value of the other individual standards and unknown sample replicates
- 8. Prepare a standard curve by plotting the average Blank-corrected 750 nm values for each BSA standard vs. its concentration in μ g/mL. Use the standard curve to determine the protein concentration for each unknown sample.

Reagents

Elution Buffer, 250 ml

In 250 ml dH₂O

- + 1.73 g NaH₂PO₄•H₂O, 50 mM (FW 137.99 g/mol)
- + 4.39 g NaCl, 300 mM (FW 58.44 g/mol)
- + 4.25 g Imidazole, 250 mM (FW 68.08 g/mol)

pH to 8.0 using NaOH store at 4°C

Allison Sutter