

PCR

Monday, December 05, 2011

1:10 PM

Sean Weise

<http://signal.salk.edu/tdnaprimers.2.html>

PCR

Materials:

Sterile dH₂O

5 U / μ l Taq (Invitrogen 18038-042)

10X PCR Buffer (comes with Invitrogen Taq)

50 mM MgCl₂ (comes with Invitrogen Taq)

2.5 mM each dNTP solution (Invitrogen 18038-042)

10 μ M Forward and Reverse primers

Agarose

TAE solution (1X and 10X)

Gel casting mold, combs, and electrophoresis apparatus

EZ-vision Loading buffer / DNA dye (Amresco, N472-Q-0.5ml)

1 Kb ladder (Amresco, K181-500 μ l)

Method:

Do the following methods on ice

1. Prepare a "Master Mix" in a 0.5 ml microfuge tube for your PCR reactions according to the following table

Reagent (conc in PCR tube)	1X (μ l)
Sterile dH ₂ O	12.5
10X PCR Buffer (1X)	2
50 mM MgCl ₂ (2.5 mM)	1
10 mM dNTPs (150 μ M)	0.3
5 U/ μ l Taq (1 U)	0.2*
Sub Total	16
10 μ M Forward Primer (1 μ M)	1
10 μ M Reverse Primer (1 μ M)	1
DNA	2
Total Volume	20

* Could go down to 0.04 - 0.1 ul if needed according to invitrogen protocol

2. Add 0.5 - 4 μ l of DNA sample to each PCR tube

0.5 μ l for purified plasmid to 4 μ l for colony PCR from liquid culture

3. Add Master Mix to each PCR tube

4. If you did not add your primers to the master mix add them now

5. If using the old Perkin Elmer thermocycler add 20 μ l of mineral oil to top

6. Place your PCR tubes in the thermocycler and run your program¹

7. While you are waiting for the PCR you can make a gel to check your results

8. In a 250 - 500 ml bottle or flask make up \approx 1 % agarose solution

For small gels add 0.5 g agarose to 50 ml 1X TAE

For large gels add 2 g agarose to 200 ml 1X TAE

9. Microwave agarose and TAE buffer for 1 minute or until boiling

10. Set up gel mold in empty tray by placing combs in mold and clamp mold tight

11. Slowly pour agarose TAE solution into gel mold, allow to cool and solidify for 30 min

12. When the agarose gel has solidified carefully and slowly remove comb and place in electrophoresis apparatus, don't worry if some of the wells are stuck closed they will open when you add TAE buffer in the next step. You can also store your gel in the refrigerator overnight. If your using a mold or comb that is in "high demand" remove the gel from the mold too. Place the gel in a plastic bag so that it does not dry out.

13. Place your gel in the electrophoresis system and fill with 1X TAE to the fill line, you can use the TAE solution that is already

in the system, no need to drain and refill.

14. When your PCR is finished, add 5 μ l of loading buffer to 20 μ l DNA in separate PCR tube
- 15 Mix 10 μ l of 1 Kb ladder with 2.5 μ l EZ-vision loading dye
16. Place 5 – 20 μ l of PCR product with loading buffer into each well.
17. After loading place cover on electrophoresis system and run gel at approximately 85 V for 30 – 45 minutes
18. Turn power supply off and remove cover. Carefully place gel in an empty tray to carry to the imager.

¹ The program you will run in the thermocycler will depend on the size of your primers and the size of the expected product, and the Taq polymerase you are using, below is an example program

Typical Program

1. 95C – 2 min ^{*1}
2. 95C – 30sec
3. 55C^{*2} – 30sec
4. 72°C – 1.5 min ^{*3}
5. Cycle to step 2 35 more times
6. 72°C – 4.5 min (3 * the extension time length)
7. Incubate at 4°C forever

^{*1} The initial melting temperature may be as long as 7.5 minutes for colony PCR to break open cells

^{*2} The annealing temperature is the most critical, a rule of thumb is to subtract 5°C from the T_m given to you by the primer 3 web program or you can check this temp by entering your primer sequence in promega's website program getting the "base stacking T_m" subtracting 2 - 3°C and using that as your annealing temperature

<http://fokker.wi.mit.edu/primer3/input.htm>

<http://www.promega.com/biomath/calc11.htm>

^{*3} The extension time should be \approx 1 kb per minute

² The % of agarose in your gel will be determined by the size of the PCR product you are trying to isolate. The smaller the product size (in base pairs) the higher the % of agarose you will need to use. In general 0.5% = 20 kbp – 2 kpb, 1% = 10kpb – 500bp, 1.5% = 5kbp – 250bp, 2% 2.5kbp – 100bp

PCR Reagents

10 X TAE Buffer 1L

- + 48.4 g Tris
- + 7.44g Na₂EDTA * 2H₂O
- + 11.42 ml glacial acetic acid
- + 100 μ l 0.5 M EDTA
- Add H₂O up to 1 L
- pH to 8.5

Note – Do not make up more than 1 L it will go bad before you can use it all

1 X TAE Buffer 1 L

- In 900 ml dH₂O
- + 100 ml 10 X TAE buffer (40 mM Tris)

100 μ M Primers \approx 1 ml

Multiply 10 * the nmol amount of primer printed on the side of the tube
This is the amount of sterile dH₂O you should add to the tube
Mix by up and down pipetting try not to vortex
store in -20°C

10 μ M Primers 500 μ l

- In 450 μ l sterile dH₂O
- + 50 μ l of 100 μ M primer
- store in -20°C