

Sandwich ELISA

2011年12月12日
15:30

Ziru Li

Sandwich ELISA protocol - need optimization

Last modified 07/25/08

1. Incubate plate with capture antibody in carbonate buffer overnight.
 - 100ul/well at 4° C; seal plate.
2. Wash 3 times with PBS.
 - 100ul/well
3. Block wells with PBS-1%M.
 - 100ul/well on shaker for 2 hours; seal plate.
4. Wash 3 times with PBS-T.
 - 100ul/well
5. Incubate wells with serial dilutions of antigen in PBS.
 - 100ul/well on shaker for 2 hours; seal plate.
6. Wash 3 times with PBS-T.
 - 100ul/well
7. Incubate wells with detection antibody in PBS-0.2%M.
 - 100ul/well on shaker for 1 hour; seal plate.
8. Wash 6 times with PBS-T.
 - 150ul/well on shaker 2min each.
9. Add substrate. [Always use a new trough for adding substrate]
 - Colorimetric (TMB): read $A_{370/655}$ on plate reader in kinetics mode right after substrate is added.
Alternative ways:
 - Read endpoint absorbance after 45 minutes incubation on shaker.
 - Add stop solution (1M H_2SO_4) at 45 minutes and read absorbance at 450nm (optimization needed)
 - Chemiluminescent (Pierce): Mix reagent 1&2 in 1:1 ratio and read on plate reader after 5 minutes incubation.

Bicarbonate buffer: 1.59g Na_2CO_3 , 2.93g $NaHCO_3$, DI H_2O to 1000mL, pH to 9.6

10xPBS: NaCl 80.0g, KCl 2g, $Na_2HPO_4 \cdot 7H_2O$ 17.3g, KH_2PO_4 2g, DI H_2O to 1000mL, pH to 7.2

PBS: 10xPBS 100mL, DI H_2O 900mL, pH to 7.2

PBS-1%BSA: 1% BSA in PBS (e.g. 2g BSA in 200mL PBS)

PBS-T: 0.05% Tween20 in PBS-1%BSA (e.g. 250uL Tween20 in 500mL PBS-1%BSA)

PBS-0.2%BSA: 0.2% BSA in PBS (e.g. 0.2g BSA in 100mL PBS)

Or

(Remake the following every week)

PBS-1%M: 1% non-fat dry milk in PBS (e.g. 2g dry milk in 200mL PBS)

PBS-T: 0.05% Tween20 in PBS-1%M (e.g. 250uL Tween20 in 500mL PBS-1%M)

PBS-0.2%M: 0.2% non-fat dry milk in PBS (e.g. 0.2g dry milk in 100mL PBS)