

# qPCR

Friday, October 28, 2011  
10:46 AM

Sean Weise

## qPCR

### Materials:

RNase Free H<sub>2</sub>O (Fisher BP5611)  
SYBR Green PCR Master Mix (ABI 4309155)  
96 Well PCR Optical plate (Eppendorf twin.tec skirted, cat. no. 951020401)  
Optical adhesive cover (Eppendorf PCR film, cat. no. 951023019)  
10 µM Forward Primer  
10 µM Reverse Primer  
cDNA Template

### Method:

1. Make up master mix on ice according to the following formula  
+ 10 µl 2X SYBR Green PCR Master Mix \* # of wells  
+ 1 µl 10 µM Forward Primer \* # of wells  
+ 1 µl 10 µM Reverse Primer \* # of wells  
+ 6 µl RNase free H<sub>2</sub>O \* # of wells
2. Pipette 18 µl of master mix into each well, keep plate cold
3. Add 2 µl cDNA Template made from 300 ng - 1 µg RNA
4. Put Optical adhesive cover on top the plate
5. Place plate in the thermocycler and run your program

### Typical Program

1. 95C – 10 min
2. 95C – 15 sec
3. 60C – 1 min
5. Cycle to step 2 40 more times
6. Melting curve
7. Incubate at 4°C forever

**Note:** When generating standards use a high fidelity DNA polymerase

### Typical Concentrations used for qPCR standards

Actin2

- For cDNA made with 500 ng RNA  
15,625,000 copies/µl  
3,125,000  
625,000  
125,000

IPP2

- For cDNA made with 300 ng RNA  
125,000 copies/µl  
25,000  
5,000  
1,000

### Math

To convert from ng/μl (obtained from spec) to copies/μl use the following formula

$$X \text{ ng}/\mu\text{l} \div 1 \cdot 10^{-9} \text{ ng}^* \div \text{molec wt of standard}^{**} \text{ g/mol} * 6.022 \times 10^{23} \text{ copies/mol} \div 2^{***} = \text{copies}/\mu\text{l}$$

\* This converts from ng/μl to g/μl

\*\* To figure out the molecular weight of your PCR product you are going to use for a standard put your sequence into this website <http://www.basic.northwestern.edu/biotools/oligocalc.html>

- be sure to select dsDNA on the drop down box on the right (default is ssDNA)

- Do Not Use Vector NTI For This!, it only gives you the molecular weight of a single stand (it's meant for use with primers)

\*\*\* We divide by 2 since RNA is single stranded i.e. 1 strand not 2

Use 100x dilutions to get PCR product into range for standards needed (i.e. 2 μl PCR product into 198 μl dH<sub>2</sub>O)